Celt ACD[®]

Vascular Closure Device

INSTRUCTION FOR USE

READ ALL INSTRUCTIONS,
WARNINGS, CAUTIONS AND SAFETY
INFORMATION CAREFULLY PRIOR TO USE

Vasorum Ltd

Trinitas House, 2012 Orchard Avenue, Citywest Campus,
Dublin 24, Ireland TEL: +353 1 4035460

CE
0050

SECTION A – DESCRIPTION

The **Celf ACD®** Vascular Closure Device is a puncture closure device which is radiopaque and visible under fluoroscopic control during all phases of deployment and implantation. It is designed to close a femoral arteriotomy following either a diagnostic or interventional catheterisation procedure. The device consists of either a 5F, 6F or 7F delivery system, with a premounted implant on its tip. Once located, the device is activated and delivers a small biocompatible stainless steel (316LVM) implant (see contraindications Section D) into the vascular wall thereby closing the arteriotomy and effecting haemostasis.

SECTION B - COMPONENTS

The **Celf ACD®** consists of the following (Fig. 1):

1. A delivery system.
2. Pre-attached implant.

SHEATH SIZE	MAX OVERALL LENGTH	MIN LUMEN DIAMETER
5F	15CM	1.7MM
6F	15CM	2.1MM
7F	15CM	2.4MM

TABLE 1

SECTION C - INTENDED USE

The **Celf ACD®** is indicated for the percutaneous closure of the common femoral arteriotomy of adult patients following diagnostic and/or interventional therapeutic intra-arterial procedures such as peripheral or coronary angiography, arterial stents, balloon angioplasty and atherectomy procedures where the arteriotomy is in the common femoral artery and introducer sheaths (either **5F**, **6F** or **7F**) have been used.

CAUTION:

The device is intended for use with the pre-existing introducer sheath. The device should only be used with sheaths that have a maximum overall length not greater than that shown in Table 1 and an internal lumen diameter which is not less than that shown in Table 1.

It is essential that the device correctly fits the intended introducer sheath.

SECTION D - CONTRAINDICATIONS

The **Celf ACD®** Vascular Closure Device is contraindicated in:

1. Patients with a femoral artery lumen diameter less than 5mm for 5F and 6F sheaths and less than 6mm for 7F sheath size.
2. Patients with a known allergy to stainless steel (316LVM) implant material.
3. Severe, acute non-cardiac systemic disease or terminal illness with a life expectancy of less than one year.
4. Evidence of systemic bacterial or cutaneous infection, including groin infection.
5. Patients suffering with definitive or potential coagulopathy or platelet count < 100,000/ μ l.
6. Patients in which introducer sheaths smaller than 5F or greater than 7F have been used.
7. Patients judged unsuitable for surgical repair or having significant peripheral vascular disease affecting the limb.
8. If puncture site is via a vascular graft.
9. If a palpable haematoma is observed during the procedure.
10. If there is any indication that the puncture has been made in the profunda femoris.
11. Patients with an active haematoma, arteriovenous fistula, or pseudoaneurysm.
12. Patients with a very superficial artery where the depth from skin to the artery surface at the access site is less than 4mm.

SECTION E - POTENTIAL COMPLICATIONS

1. Bruising, oozing or bleeding at the puncture site.
2. Haematoma or Ecchymosis.
3. Pain, discomfort or transitory local irritation and inflammation at the puncture site.
4. Access-site related nerve injury, vascular spasm, local pulse deficits, ischaemia or access-site wound dehiscence.
5. Infection at the procedure site.
6. Arterial or deep vein thrombosis, pseudoaneurysm or arteriovenous fistula.
7. Failure of the device to deploy correctly in the artery.
8. Vaso-vagal response is a common and expected response to manipulation of the blood vessels and may occur during any invasive vascular procedure.

Other events that could possibly occur, but are considered extremely unlikely include:

9. Occlusive intraluminal thrombus and/or emboli formation at the **Celf ACD®** Vascular Closure Device implant site.
10. Sinus formation.
11. Partial or complete occlusion of the arterial lumen.
12. Damage of the arterial wall including stenosis/narrowing.
13. Swelling of the treated limb.
14. Corrective intervention due to any of the above complications such as surgery and/or transfusion.
15. Migration of the implant.
16. Pseudoaneurysm.

SECTION F - WARNINGS

1. Do not use if package is damaged or any portion of the package has been previously opened.
2. Do not use if the label is missing or not indicating that the device is sterile.
3. Do not use any of the components if they appear damaged or defective in any way.
4. The **Celf ACD**[®] Vascular Closure Device represents a foreign body and therefore it is not recommended to use the device in situations where bacterial contamination of the surrounding tissue may have occurred.
5. The **Celf ACD**[®] Vascular Closure Device is intended for single use and cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse. Attempts to clean or sterilize these devices may result in bio-incompatibility, infection or produce failure risks to the patient.

SECTION G - PRECAUTIONS

1. The **Celf ACD**[®] Vascular Closure Device is intended for use by healthcare professionals and personnel possessing adequate instruction in the use of the device and trained in intravascular catheter procedures.
2. A single wall puncture technique must be used to ensure that puncture of the posterior wall of the artery does not occur.
3. Aseptic sterile techniques should be observed at all times when using the **Celf ACD**[®] Vascular Closure Device.
4. Exercise care during device handling to reduce the possibility of accidental damage to the **Celf ACD**[®] Vascular Closure Device.
5. In the event of persistent bleeding from the femoral access site after use of the **Celf ACD**[®] Vascular Closure Device, use standard manual or mechanical compression techniques.
6. Disposal of contaminated devices, components, and packaging materials should follow universal precautions for bio-hazardous waste.

SECTION H – PROCEDURE

Indicator Arrows:

Numbered indicator arrows are printed onto the body of the handle to ensure that the device is deployed correctly (Fig. 2).

Arrow 1 is printed onto the handle indicating that it is the first step in the deployment sequence. The handle should be rotated in a clockwise direction (indicated by the arrow direction) to deploy the distal wings (Step 3).

Arrow 2 is also printed onto the handle indicating that the handle should be turned in an anti-clockwise direction to deploy the proximal wings (Step 5).

Arrow 3 is printed onto both sides of the ejector lever indicating the direction in which the lever should be pressed to eject the implant (Step 7). Additional information such as the French size and the trade name for the device (**Celf ACD**[®]) are also printed onto the device.

Step 1:

- a. Refer to the pre-procedure fluoroscopic image A below to confirm the site of puncture by the sheath into the artery.
- b. Under strict sterile conditions and using an aseptic technique, remove the **Celf ACD**[®] Vascular Closure Device from its packaging.
- c. Insert the tip of the **Celf ACD**[®] Vascular Closure Device through the valve in the introducer sheath and advance the flexible guide tube into the lumen of the sheath (Fig. 3).

Step 2:

- a. Continue to advance the **Celf ACD**[®] Vascular Closure Device through the introducer sheath and into the femoral artery until the latches extending from the front of the handle are touching the valve of the sheath. Fluoroscopic control is used to confirm position of the device within the arterial lumen.
- b. Holding the device in a fixed position, pull the valve of the introducer sheath toward the handle of the delivery system and into the latches extending from its front end. Ensure that the valve of the introducer sheath is securely attached to the device (Fig. 4).

IMPORTANT: Do not continue to deploy the implant unless the Celf ACD[®] Vascular Closure Device is properly engaged with the introducer sheath.

Step 3:

- a. Holding the device and the introducer sheath, withdraw both from the patient until the distal tip of the sheath is close to the puncture hole but remains within the arterial lumen. Fluoroscopic control may be used for guidance.
- b. Prior to operating the device, remove the red locking tab by pressing down on the locking tab lever (Fig. 5).
- c. Holding the **Celf ACD**[®] Vascular Closure Device with one hand, grip the end of the handle with the other hand and turn it **GENTLY** in a **clockwise direction**, as indicated by **Arrow 1**, until it pops back to a stop position (Fig. 6). The opening of the distal wings of the implant may be observed under fluoroscopy.

IMPORTANT: Check to ensure that the handle has moved back leaving a gap, as illustrated below (Fig. 7), before proceeding to the next step.

Step 4:

- a. Holding the **Celf ACD**[®] Vascular Closure Device with one hand and the valve of the introducer sheath with the other hand, withdraw it from the patient until resistance is felt (Fig. 8).

- b. If fluoroscopic image taken, the position of the distal wings at the puncture wall should correlate with the puncture sheath entry site as seen on the pre deployment fluoroscopic image.

IMPORTANT: Ensure that the introducer sheath remains attached to the handle during withdrawal (Fig. 8).

Step 5:

- a. To deploy the proximal wings, firstly check to ensure that the valve of the sheath remains attached to the device. **Hold the device as close to vertical as possible (Fig. 9).**
- b. While maintaining upward tension on the device and the vertical orientation, turn the handle **GENTLY** in a **counter-clockwise** direction, as indicated by **Arrow 2**, until it pops back to the second stop position (Fig. 9).
- c. Fluoroscopic control may be used to observe the proximal wings opening on the external surface of the artery to form a "sandwich" with arterial puncture edges between the wings.

IMPORTANT: Check to ensure that the handle has moved back leaving an increased gap before proceeding to the next step (Fig. 9).

Step 6:

- a. Carry out a sandwich "Push-Pull" test by gently pushing the device backward and forward. If deployment is correct, it will not be possible to advance the deployed implant into the artery.
- b. If fluoroscopic image taken, the implant will be observed to remain in a fixed position.

Step 7:

- a. It is now safe to eject the implant from the delivery system by pressing down on the ejector lever, as indicated by **Arrow 3**, until it abuts the handle (Fig. 10).
- b. If fluoroscopic image taken, the separation of the delivery device from the implant may be observed.

Step 8:

- a. Remove the device and introducer sheath from the patient and dispose of both appropriately (Fig. 11).
- b. If required a final fluoroscopic picture will show the implant stabilised in the arterial wall.

SECTION I - CAUTION

- a. Each stage of the deployment has been designed to proceed in a smooth and orderly manner. If resistance is felt at any stage, this means that the previous part of the cycle has not been fully completed.
- b. In the event that post-firing excessive bleeding occurs with the removal of the **Celf ACD®** Vascular Closure Device, apply either manual or mechanical pressure as would normally be used when an arterial puncture closure device is not in use.

MRI Safety Information:

Non-Clinical testing demonstrated that the **Celf ACD®** Vascular Closure Device is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less.
- Maximum spatial gradient magnetic field of 1,500-gauss/cm (15-T/m) (extrapolated) or less.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the First Level Controlled Operating Mode of operation for the MR System.

Under the scan conditions defined for the **Celf ACD®** is expected to produce a maximum temperature rise of 2.1°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artefact caused by the **Celf ACD®** extends approximately 10-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

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